




# THE FEATURES

**Easy to use**  
2 min hands-on time

**Fast and labor-saving**  
result in ~100 min for 19 pathogens

## Interconnected Tube




**Safe**  
disposable sealed cartridge

**Ambient transport**  
cartridge pre-filled with glassified reagent,  
no need to redissolve it before use

**Point-of-care testing**  
independent modules  
run any test any time with 24/7 random access

## Helicobacter pylori 23S rRNA Gene mutation Detection Kit

Target genes	<ul style="list-style-type: none"> <li><i>Helicobacter Pylori</i> (UreA)</li> <li>Clarithromycin resistance genes: 23S rRNA genes A2142G, A2142C, A2143G</li> <li>Quinololon resistance genes: gyrA genes A260T, C261A, T261G, G271A, G271T, A272G</li> </ul>
LoD	1000 copies/ml
Time to result	90 min
Sample material	Biopsy



**Ustar Biotechnologies (Hangzhou) Ltd.**