

EasyNAT



Introducing the Ustar EasyNAT® System

CE-IVDR

Experience a revolutionary platform that brings laboratory-quality nucleic acid testing to the point of care. The Ustar EasyNAT® system features a growing CE-IVDR menu for unparalleled convenience and accuracy.

The Power of EasyNAT®

Speed | Accuracy | Flexibility

The Ustar EasyNAT® system delivers unmatched results at the point of care with remarkable speed, precision, and versatility.



• CE-IVDR Test Menu •

Respiratory Health:

- EasyNAT Rapid Flu Assay
- EasyNAT Rapid RSV Assay
- EasyNAT MP Assay
- EasyNAT BP Assay

Critical Infectious Disease:

- EasyNAT MTC Assay

Women's & Sexual Health:

- EasyNAT GBS Assay
- EasyNAT CT/NG Assay
- EasyNAT MG Assay
- EasyNAT HSV 1&2 Assay
- EasyNAT UU Assay

Designed for Point of Care Excellence



Compact Design: Reduced footprint with built-in barcode scanner and full-size touchscreen.



Ultimate Flexibility: On-demand and random access with mix-and-match test capabilities.



User-Friendly: No specialized training or daily maintenance required.



Effortless Printing: Simple, remote printing capabilities.



Network Ready: Flexible wired and wireless network connection capabilities.


Easy as 1,2,3

- With less than one minute of hands-on time

1 Load Sample



2 Insert Cartridge



3 Report Result



Specifications

CATALOG NUMBER	DESCRIPTION	INSTRUMENT DIMENSIONS
UC0204	4-module system	300 mm (L) x 470 mm (W) x 390 mm (H)
UC0208	8-module system	300 mm (L) x 470 mm (W) x 390 mm (H)
UC0216	16-module system	570 mm (L) x 470 mm (W) x 390 mm (H)

Simplify Your Workflow

- **Space-Efficient:** Compact instrument runs up to 16 tests simultaneously, saving valuable bench space.
- **User-Friendly Interface:** Intuitive software and technology minimize training needs.
- **Fully Automated:** Sample extraction, nucleic acid amplification, and detection are all automated.
- **Safe and Clean:** Self-contained cartridge reduces waste and contamination risks.
- **Built-In Quality Control:** Integrated quality control in every cartridge.
- **LIS Integration:** Capabilities for interfacing with Laboratory Information Systems (LIS).

Discover how the Ustar EasyNAT® System can transform your point of care testing.

EasyNAT Tests

Infections	Tests	Catalog Number	Number of Tests	Certification
Respiratory	EasyNAT SARS-CoV-2 Assay	U20223	20 tests/kit	CE-IVD
	EasyNAT Rapid Flu Assay	U20206	20 tests/kit	CE-IVDR
	EasyNAT Rapid RSV Assay	U20208	20 tests/kit	CE-IVDR
	EasyNAT MP Assay	U20201	20 tests/kit	CE-IVDR
	EasyNAT BP Assay	U20209	20 tests/kit	CE-IVDR
	EasyNAT AdV/HRV Assay	U202059-20	20 tests/kit	CE-IVD
	EasyNAT SP Assay	U202005-20	20 tests/kit	CE-IVD
	EasyNAT LP Assay	U202006-20	20 tests/kit	CE-IVD
Women's Health, & Sexual Health	EasyNAT GBS Assay	U20511	20 tests/kit	CE-IVDR
	EasyNAT CT/NG Assay	U205001-20	20 tests/kit	CE-IVDR
	EasyNAT MG Assay	U20502	20 tests/kit	CE-IVDR
	EasyNAT HSV 1&2 Assay	U20504	20 tests/kit	CE-IVDR
	EasyNAT UU Assay	U20503	20 tests/kit	CE-IVDR
	EasyNAT UU/MH Assay	U205005-20	20 tests/kit	CE-IVD
	EasyNAT TV Assay	U205009-20	20 tests/kit	CE-IVD
	EasyNAT TP Assay	U205004-20	20 tests/kit	CE-IVD
	EasyNAT Human Parvovirus B19 Assay	U205002-20	20 tests/kit	CE-IVD
	EasyNAT HPV 6/11 Assay	U205012-20	20 tests/kit	CE-IVD
	EasyNAT HPV 16/18 Assay	U205011-20	20 tests/kit	CE-IVD
TB & Emerging Infectious Diseases	EasyNAT MTC Assay	U20103	20 tests/kit	CE-IVDR
	EasyNAT MTC/NTM Assay	U201002-20	20 tests/kit	CE-IVD
	EasyNAT Monkeypox Virus Assay	U202028-20	20 tests/kit	CE-IVD
	EasyNAT Malaria Assay	U203010-20	20 tests/kit	CE-IVD

Visit our website to discover how the Ustar EasyNAT® System can transform your point of care testing.

Add: 6/F, Building 2, 611 Dongguan Road, Binjiang District, Hangzhou, Zhejiang, 310053 PR China
 Tel: +86-571-88939358
 Email: sales@bioustar.com
 Web: en.bioustar.com

Ustar Biotechnologies (Hangzhou) Ltd.



Facebook

LinkedIn

YouTube